

Hemoglobin alpha antibody [EG1-8]

Cat. No. GTX04752

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P, FCM, IP
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
IHC-P	1:50-1:200
FCM	1:50
IP	1:50

Not tested in other applications.

Calculated MW 15 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 150 mM NaCl, 50% glycerol, 0.4-0.5 mg/mL BSA
Preservative	0.02% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthesized peptide derived from human Hemoglobin subunit alpha.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04752 IHC-P Image

IHC-P analysis of human spleen tissue using GTX04752 Hemoglobin alpha antibody [EGI-8].


GTX04752 IHC-P Image

IHC-P analysis of mouse embryo tissue using GTX04752 Hemoglobin alpha antibody [EGI-8].
Antigen retrieval : Heat mediated Tris-EDTA buffer (pH 8.0)


GTX04752 WB Image

WB analysis of K562 whole cell lysate using GTX04752 Hemoglobin alpha antibody [EGI-8].



For full product information, images and publications, please visit our [website](https://www.genetex.com).