

EZH2 antibody [IDF-5]

Cat. No. GTX04778

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:50-1:200
IHC-P	1:50-1:200
FCM	1:50

Not tested in other applications.

Calculated MW 85 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 150 mM NaCl, 50% glycerol, 0.4-0.5 mg/mL BSA
Preservative	0.02% sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthesized peptide derived from human KMT6 / EZH2.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04778 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX04778 EZH2 antibody [IDF-5].


GTX04778 IHC-P Image

IHC-P analysis of human tonsil tissue using GTX04778 EZH2 antibody [IDF-5].


GTX04778 WB Image

WB analysis of HEK293 whole cell lysate using GTX04778 EZH2 antibody [IDF-5].



For full product information, images and publications, please visit our [website](https://www.genetex.com).