

PAX2 antibody

Cat. No. GTX04784

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.5 µg/mL
IHC-P	0.5-1 µg/mL

Note : By heat

Not tested in other applications.

Calculated MW 45 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	5 mg BSA, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄
Preservative	0.05 mg sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Pax2 protein.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx04784 WB Image


WB analysis of various samples using GTx04784 PAX2 antibody.

Lane 1 : Rat brain tissue lysate

Lane 2 : Rat lung tissue

Loading : 50 µg

Dilution : 0.5 µg/mL




GTx04784 IHC-P Image

IHC-P analysis of rat lymphaden tissue using GTx04784 PAX2 antibody.

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)

Dilution : 1 µg/mL




GTx04784 IHC-P Image

IHC-P analysis of mouse lymphaden tissue using GTx04784 PAX2 antibody.

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)

Dilution : 1 µg/mL



GTx04784 IHC-P Image

IHC-P analysis of human tonsil tissue using GTx04784 PAX2 antibody.

Antigen retrieval : Heat mediated EDTA buffer (pH 8.0)

Dilution : 1 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).