

PAX2 antibody

Cat. No. GTX04786

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human, Mouse, Rat, Sheep

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
ICC/IF	Assay dependent
IHC-P	Assay dependent
FCM	1:100

Note : Permeabilization : 90% ice-cold methanol

Not tested in other applications.

Calculated MW 45 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.01M TBS, 1% BSA, 50% Glycerol
Preservative	0.02% Proclin300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human PAX2.
Purification	Protein A purified
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTx04786 IHC-P Image

IHC-P analysis of rat lung tissue using GTx04786 PAX2 antibody.

Dilution : 1:200



GTx04786 FCM Image

FACS analysis of 4% PFA-fixed HepG2 cells using GTx04786 PAX2 antibody.


Green : Primary antibody

Blue : Isotype control

Black : Unlabeled samples

Permeabilization : 90% ice-cold methanol

Dilution : 1:100



GTx04786 IHC-P Image

IHC-P analysis of mouse pancreas tissue using GTx04786 PAX2 antibody.

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).