

## Streptococcus Group A antibody [C01756R]

Cat. No. GTX04789

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	Lateral Flow, Sandwich ELISA
Reactivity	Streptococcus Group A

Package  
500 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Lateral Flow	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture : GTX04789, Detection : GTX04790 or Capture : GTX04790, Detection : GTX04789**

Not tested in other applications.

**Product Note** No cross reactivity with Streptococci Groups B, C, D, F, and G and other bacteria.

## Properties

Form	Liquid
Buffer	0.01 M PBS
Preservative	0.1% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	2.53 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Streptococcus Group A polysaccharide.
Purification	Protein A purified
Purity	>90% by SDS-PAGE
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).