

Streptococcus Group A antibody

Cat. No. GTX04790

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	Sandwich ELISA
Reactivity	Streptococcus Group A

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Sandwich ELISA	Assay dependent

Note : Capture : GTX04789, Detection : GTX04790 or Capture : GTX04790, Detection : GTX04789

Not tested in other applications.

Product Note Type specific carbohydrate of group A Streptococcus. Does not react with other Strep groups.

Properties

Form	Liquid
Buffer	0.01 M PBS
Preservative	0.1% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	5.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Streptococci, group A
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).