

Rabbit Anti-Mouse IgG antibody, pre-adsorbed (AP)

Cat. No. GTX04807

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA, IHC
Reactivity	Mouse

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
ELISA	1:9800
IHC	1:100-1:500

Not tested in other applications.

Product Note No reactivity was observed against Human serum proteins.

Properties

Form	Liquid
Buffer	50mM Tris-HCl, 150mM NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂ , 50% Glycerol, 1%BSA
Preservative	0.01% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse IgG whole molecule
Purification	Purified by antigen-affinity chromatography followed by solid phase adsorption to remove unwanted reactivities and extensive dialysis against the buffer From serum
Conjugation	Alkaline Phosphatase (AP)

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx04807 WB Image

WB analysis of Mouse IgG protein using GTx04807 Rabbit Anti-Mouse IgG antibody, pre-adsorbed (AP).

Loading : 0.1 µg

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).