

Rabbit Anti-Human IgG (Fc) antibody (HRP)

Cat. No. GTX04809

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA, IHC
Reactivity	Human

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:5000
ELISA	1:10000-1:50000
IHC	1:500-1:2500

Not tested in other applications.

Properties

Form	Liquid
Buffer	10mM Sodium Phosphate, 150mM NaCl, 1% BSA
Preservative	0.01% Gentamicin Sulfate. Do NOT add Sodium Azide!
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human IgG F(c) fragment
Purification	Purified by chromatography From serum
Conjugation	Horseradish peroxidase(HRP)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04809 ELISA Image

ELISA analysis of purified Human IgG (Fc) protein using serial diluted GTX04809 Rabbit Anti-Human IgG (Fc) antibody (HRP). The starting dilution of antibody was 5µg/ml.

Blue line: 1µg Human IgG (Fc) protein



For full product information, images and publications, please visit our [website](https://www.genetex.com).