

Plasmodium falciparum pLDH antibody [3560]

Cat. No. GTX04821

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, ELISA, Sandwich ELISA
Reactivity	Plasmodium falciparum

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:10-1:50
ELISA	1:20-1:200
Sandwich ELISA	Assay dependent

Note : Acetone fixation is recommended prior to staining
Capture : GTX04822, Detection : GTX04823.

Not tested in other applications.

Product Note

Specific for the parasitic lactate dehydrogenase (pLDH) of several Plasmodium species: Plasmodium falciparum, Plasmodium malariae and Plasmodium vivax.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	full-length recombinant LDH protein
Purification	Protein A purified
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 2

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).