

NKX6-1 antibody [NKX61/2561]

Cat. No. GTX04829

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2c
Applications	IHC-P, Protein Array
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1- 2 µg/mL
Protein Array	Assay dependent

Note : Recommend using heat mediated antigen retrieval in 10mM Tris with 1mM EDTA (pH 9.0)

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA (Please contact us for PBS only format)
Preservative	0.05% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full-length human NKX6.1 protein
Purification	Protein A/G purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04829 Protein Array Image

Analysis of Protein Array containing more than 19,000 full-length human proteins using NKX6-1 antibody [NKX61/2561]. Z- and S- Score: The Z-score represents the strength of a signal that a monoclonal antibody produces when binding to a particular protein on the HuProt™ array. Z-scores are described in units of standard deviations (SD's) above the mean value of all signals generated on that array. If targets on HuProt™ are arranged in descending order of the Z-score, the S-score is the difference (also in units of SD's) between the Z-score. S-score therefore represents the relative target specificity of a Monoclonal Antibody to its intended target. A Monoclonal Antibody is considered to specific to its intended target if the Monoclonal Antibody has an S-score of at least 2.5. For example, if a Monoclonal Antibody binds to protein X with a Z-score of 43 and to protein Y with a Z-score of 14, then the S-score for the binding of that Monoclonal Antibody to protein X is equal to 29.


GTx04829 IHC-P Image

IHC-P analysis of human pancreas tissue using GTx04829 NKX6-1 antibody [NKX61/2561].


GTx04829 IHC-P Image

IHC-P analysis of human pancreas tissue using GTx04829 NKX6-1 antibody [NKX61/2561].



For full product information, images and publications, please visit our [website](https://www.genetex.com).