

VAT1L antibody

Cat. No. GTX04840

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
IHC-P	1:100-1:500

Not tested in other applications.

Calculated MW 46 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	10mM TBS, 1% albumin, 50% Glycerol
Preservative	0.03% Proclin300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human VAT1L (321-419/419 aa)
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX04840 WB Image

WB analysis of various samples using GTX04840 VAT1L antibody.


Lane 1 : Mouse cerebrum tissue lysate

Lane 2 : Mouse cerebellum tissue lysate

Lane 3 : Rat cerebrum tissue lysate

Lane 4 : Rat cerebellum tissue lysate


Dilution : 1:1000



GTX04840 IHC-P Image

IHC-P analysis of rat brain tissue using GTX04840 VAT1L antibody.

Dilution : 1:200



GTX04840 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX04840 VAT1L antibody.

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).