

LAPTM4A antibody

Cat. No. GTX04850

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:300-1:5000

Not tested in other applications.

Calculated MW 27 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	10mM TBS, 1% BSA, 50% Glycerol
Preservative	0.02% Proclin300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human LAPTM4A (101-200/233 aa)
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04850 WB Image

WB analysis of mouse lung tissue lysate using GTX04850 LAPT M4A antibody.

Dilution : 1:1000



GTX04850 WB Image

WB analysis of various samples using GTX04850 LAPT M4A antibody.

Lane 1 : Mouse uterus tissue lysate

Lane 2 : Mouse ovary tissue lysate

Dilution : 1:1000


GTX04850 WB Image

WB analysis of mouse lung tissue lysate using GTX04850 LAPT M4A antibody.

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).