

Enterovirus VP3 antibody [C01700M]

Cat. No. GTX04852

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ICC/IF, ELISA
Reactivity	Enterovirus

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Product Note

Reacts with 9 Rhinoviruses, 6 Coxsackie B viruses, 6 ECHO viruses, 4 Coxsackie A viruses, Enterovirus type 71 and Poliovirus. It does not react with Influenza A virus, Influenza B virus, Respiratory Syncytial Virus, Adenovirus, Parainfluenza virus 1, Parainfluenza virus 2, Parainfluenza virus 3, Chlamydia pneumoniae, Mycoplasma pneumoniae, and Hepatitis A virus.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	2.7 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant VP3 from Rhinovirus type 14.
Purification	Protein A purified
Purity	> 90%
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).