

CD8 alpha antibody

Cat. No. GTX04862

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Rat, Chicken

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:300-1:5000
IHC-P	1:200-1:400

Not tested in other applications.

Calculated MW 26 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	10mM TBS, 1% BSA, 50% Glycerol
Preservative	0.02% Proclin300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from rat CD8 alpha(51-150/236 aa)
Purification	Purified by Protein A
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX04862 IHC-P Image

IHC-P analysis of rat spleen tissue using GTX04862 CD8 alpha antibody.

Antigen retrieval : Boiling sodium citrate buffer (pH6.0) for 15min

Dilution : 1:200



GTX04862 IHC-P Image


IHC-P analysis of rat spleen tissue using GTX04862 CD8 alpha antibody.

Orange : Primary antibody

Blue : DAPI

Antigen retrieval : Boiling sodium citrate buffer (pH6.0) for 15min

Dilution : 1:200




GTX04862 IHC-P Image

IHC-P analysis of rat colon tissue using GTX04862 CD8 alpha antibody.

Orange : Primary antibody

Blue : DAPI

Dilution : 1:200



GTX04862 WB Image

WB analysis of various rat tissue samples using GTX04862 CD8 alpha antibody.

Dilution : 1:300



For full product information, images and publications, please visit our [website](https://www.genetex.com).