

NKX2.8 antibody

Cat. No. GTX04864

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent

Not tested in other applications.

Calculated MW 26 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide directed towards the C terminal region of human NKX2-8
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX04864 WB Image

WB analysis of human heart tissue lysate using GTX04864 NKX2.8 antibody.
Dilution : 0.2-1 µg/mL




GTX04864 WB Image

WB analysis of human fetal lung tissue lysate using GTX04864 NKX2.8 antibody.
Dilution : 1 µg/mL



GTX04864 WB Image

WB analysis of human fetal heart tissue lysate using GTX04864 NKX2.8 antibody.
Dilution : 1 µg/mL



GTX04864 WB Image

WB analysis of mouse heart tissue lysate using GTX04864 NKX2.8 antibody.
Dilution : 1 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).