

TIM-3 antibody [1827CT151.74.68]

Cat. No. GTX04870

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, FCM
Reactivity	Human

Package
200 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:4000
FCM	1:25

Note : Permeabilization : 90% methanol for 10 min

Not tested in other applications.

Calculated MW 33 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein of human HAVCR2.
Purification	Purified by Protein G followed by dialysis against PBS
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx04870 WB Image

WB analysis of various samples using GTx04870 TIM-3 antibody [1827CT151.74.68].


Lane 1 : HepG2 whole cell lysate

Lane 2 : Daudi whole cell lysate

Lane 3 : Human kidney lysate Lysate

Loading : 20 µg

Dilution : 1:4000



GTx04870 FCM Image

FACS analysis of 2% PFA fixed U2OS cells using GTx04870 TIM-3 antibody [1827CT151.74.68].

Green line : Primary antibody

Blue line : Isotype control

Dilution : 1:25



For full product information, images and publications, please visit our [website](https://www.genetex.com).