

KLRC1 antibody

Cat. No. GTX04876

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human

Package
200 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
ICC/IF	1:25
IHC-P	1:50-1:100
FCM	1:25

Note : Permeabilization step is needed.

Not tested in other applications.

Calculated MW 26 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 180-206 amino acids from the C-terminal region of human KLRC1.
Purification	Purified by Protein A followed by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTx04876 FCM Image

FACS analysis of 2% PFA fixed Jurkat cells using GTx04876 KLRC1 antibody.


Green line : Primary antibody

Blue line : Isotype control

Dilution : 1:25


GTx04876 IHC-P Image

IHC-P analysis of human kidney carcinoma tissue using GTx04876 KLRC1 antibody.


GTx04876 WB Image

WB analysis of Jurkat whole cell lysate using GTx04876 KLRC1 antibody.

Loading : 20 µg

Dilution : 1:2000



For full product information, images and publications, please visit our [website](https://www.genetex.com).