

NR0B2 antibody [OTI7B5]

Cat. No. GTX04884

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 28 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length human recombinant protein of human NR0B2 (NP_068804) produced in E.coli.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx04884 IHC-P Image

IHC-P analysis of human liver carcinoma tissue using GTx04884 NR0B2 antibody [OTI7B5].

Antigen retrieval : Heat mediated Tris-EDTA, pH8.0


Dilution: 1:150


GTx04884 IHC-P Image

IHC-P analysis of human liver tissue using GTx04884 NR0B2 antibody [OTI7B5].

Antigen retrieval : Heat mediated Tris-EDTA, pH8.0


Dilution: 1:150


GTx04884 IHC-P Image

IHC-P analysis of human breast adenocarcinoma tissue using GTx04884 NR0B2 antibody [OTI7B5].

Antigen retrieval : Heat mediated Tris-EDTA, pH8.0


Dilution: 1:150


GTx04884 IHC-P Image

IHC-P analysis of human endometrium tissue using GTx04884 NR0B2 antibody [OTI7B5].

Antigen retrieval : Heat mediated Tris-EDTA, pH8.0

Dilution: 1:150


GTx04884 IHC-P Image

IHC-P analysis of human bladder tissue using GTx04884 NR0B2 antibody [OTI7B5].

Antigen retrieval : Heat mediated Tris-EDTA, pH8.0

Dilution: 1:150



For full product information, images and publications, please visit our [website](https://www.genetex.com).