

FOXD1 antibody

Cat. No. GTX04885

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
ICC/IF	1:25
IHC-P	1:10-1:50

Not tested in other applications.

Calculated MW 46 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human FOXD1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04885 ICC/IF Image


ICC/IF analysis of 4%PFA-fixed HepG2 cells using GTX04885 FOXD1 antibody.

Green : Primary antibody

Red : Actin

Permeabilization : 0.1% Triton X-100


Dilution : 1:25


GTX04885 WB Image

WB analysis of various samples using GTX04885 FOXD1 antibody.

Loading : 20 µg

Dilution: 1:2000


GTX04885 IHC-P Image

IHC-P analysis of human kidney tissue using GTX04885 FOXD1 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).