

Mouse IgG protein

Cat. No. GTX04893-pro

Species Mouse

Package
1 mg

Applications

Observed MW (kDa) 27|55 kDa.

Properties

Form Lyophilized powder

Buffer Reconstitute with PBS to 0.1-1mg/ml. Lyophilized from PBS, 5%Trehalose, 1mM DTT. (Please contact us for endotoxin-free format)

Preservative Proclin300

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Region/Sequence Native protein from normal mouse serum

Expression System mouse

Purification Purified by salt co-precipitation and ionic-Exchange chromatography
From serum

Purity >94% as determined by SDS-PAGE

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx04893-pro Image

SDS-PAGE of GTx04893-pro Mouse IgG protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).