

AHNAK2 antibody

Cat. No. GTX04899

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-2000
IHC-P	1:50-1:200

Not tested in other applications.

Calculated MW 617 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	10mM TBS, 1% BSA, 50% Glycerol
Preservative	0.02% Proclin300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human AHNAK2(1201-1300aa)
Purification	Purified by Protein A
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTx04899 WB Image

WB analysis of mouse cerebrum tissue lysate using GTx04899 AHNAK2 antibody.

Dilution : 1:1000



GTx04899 IHC-P Image

IHC-P analysis of rat brain tissue using GTx04899 AHNAK2 antibody.

Antigen retrieval : High pressure antigen retrieval with sodium citrate, pH6.0

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).