

p75 NGF Receptor / CD271 antibody

Cat. No. GTX04906

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:50-1:200

Not tested in other applications.

Calculated MW 45 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 273-427 of human p75NTR/NGFR (NP_002498.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX04906 ICC/IF Image

ICC/IF analysis of U-251 MG cells using GTX04906 p75 NGF Receptor / CD271 antibody.

Orange : Primary antibody

Blue : DAPI


Dilution : 1:100


GTX04906 WB Image

WB analysis of rat brain tissue lysate using GTX04906 p75 NGF Receptor / CD271 antibody.

Loading : 25µg

Dilution : 1:1000


GTX04906 WB Image

WB analysis of various samples using GTX04906 p75 NGF Receptor / CD271 antibody.

Loading : 25µg

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).