

# TRP1 antibody

**Cat. No. GTX04910**

|                  |            |
|------------------|------------|
| <b>Host</b>      | Rabbit     |
| <b>Clonality</b> | Polyclonal |
| <b>Isotype</b>   | IgG        |

**Package**  
100 µl

## Applications

**Calculated MW** 61 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                                                                                                                                      |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                                                                                                                               |
| <b>Buffer</b>        | PBS, 150mM NaCl, 50% Glycerol                                                                                                                                                                                                                                                                                                        |
| <b>Preservative</b>  | 0.02% Sodium Azide                                                                                                                                                                                                                                                                                                                   |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                           |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                             |
| <b>Immunogen</b>     | A synthesized peptide derived from human TYRP1(Accession P17643), corresponding to amino acid residues F244-E294.                                                                                                                                                                                                                    |
| <b>Purification</b>  | From serum                                                                                                                                                                                                                                                                                                                           |
| <b>Purity</b>        | Purified by antigen-affinity chromatography                                                                                                                                                                                                                                                                                          |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                                                                                                                                         |
| <b>Note</b>          | For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.<br><br>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product. |



For full product information, images and publications, please visit our [website](#).