

TIF1 alpha antibody

Cat. No. GTX04914

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, FCM, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.25 µg/ml
FCM	1-3 µg/1x10 ⁶ cells
ELISA	0.1-0.5 µg/ml

Note : Permeabilization step is needed.

Not tested in other applications.

Calculated MW 117 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄
Preservative	0.01 mg Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	E.coli-derived human TRIM24 recombinant protein (Position: A706-D964).
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04914 FCM Image


FACS analysis of 4% PFA-fixed A431 cells using GTX04914 TIF1 alpha antibody.

Blue line : Primary antibody

Green line : Isotype Control

Red line : Blank control

Dilution : 1 μ g/1x10⁶ cells


GTX04914 WB Image

WB analysis of various samples using GTX04914 TIF1 alpha antibody.

Loading : 50 μ g

Dilution : 0.25 μ g/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).