

Apolipoprotein M antibody [ARC1210]

Cat. No. GTX04920

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:50-1:200
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 21 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ApoM (O95445).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx04920 IHC-P Image

IHC-P analysis of human liver cancer tissue using GTx04920 Apolipoprotein M antibody [ARC1210].

Antigen retrieval : High pressure antigen retrieval with 0.01 M citrate buffer (pH 6.0)

Dilution : 1:100




GTx04920 WB Image

WB analysis of various samples using GTx04920 Apolipoprotein M antibody [ARC1210].

Loading : 25µg

Dilution :1:1000



GTx04920 ICC/IF Image

ICC/IF analysis of HepG2 cells using GTx04920 Apolipoprotein M antibody [ARC1210].

Red : Primary antibody

Green : α-Tubulin

Blue : DAPI

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).