

SARS-CoV-2 (COVID-19) 3CL protease protein

Cat. No. GTX04930-pro

Species SARS Coronavirus 2**Package**
100 µg

Applications

Observed MW (kDa) 38 kDa.

Properties


Form Liquid**Buffer** 20mM HEPES, 120mM NaCl, 10% Glycerol, 2mM TCEP**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** Batch dependent (Please refer to the vial label for the specific concentration.)**Region/Sequence** N-terminus His and Avi tagged; Ser1-Gln306 (YP_009725301.1)**Expression System** E. coli**Purity** >95% by SDS-PAGE**Conjugation** Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX04930-pro Image

SDS-PAGE of GTX04930-pro SARS-CoV-2 (COVID-19) 3CL protease protein.

For full product information, images and publications, please visit our [website](https://www.genetex.com).