

Integrin alpha 6 antibody [DHG-9]

Cat. No. GTX04931

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
ICC/IF	1:50-1:200
IHC-P	1:50-1:100

Not tested in other applications.

Calculated MW 127 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 150mM NaCl, 0.04-0.05% BSA, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthesized peptide derived from human Integrin alpha 6
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX04931 ICC/IF Image


ICC/IF analysis of PFA-fixed HeLa cells using GTX04931 Integrin alpha 6 antibody [DHG-9].

Green : Primary antibody

Blue : DAPI

Antigen retrieval : IHC enzyme antigen retrieval reagent


Dilution : 1:50



GTX04931 WB Image

WB analysis of SW480 whole cell lysate using GTX04931 Integrin alpha 6 antibody [DHG-9].

Loading : 50 µg



GTX04931 WB Image

WB analysis of various samples using GTX04931 Integrin alpha 6 antibody [DHG-9].

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).