

FOXA3 antibody

Cat. No. GTX04941

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.5µg/ml

Not tested in other applications.

Calculated MW 37 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 5mg BSA
Preservative	0.05mg Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human FOXA3.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX04941 WB Image

WB analysis of various samples using GTX04941 FOXA3 antibody.

Lane 1 : Rat liver tissue lysate (50 µg)

Lane 2 : Rat pancreas tissue lysate (50 µg)

Lane 3 : Mouse liver tissue lysate (50 µg)

Lane 4 : HeLa whole cell lysate (40 µg)

Dilution : 0.5 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).