

Twist1/2 antibody

Cat. No. GTX04942

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.25-0.5µg/ml
IHC-P	2-5µg/ml
FCM	1-3 µg/1x10 ⁶ cells

Note : Permeabilization step is needed.

Not tested in other applications.

Properties

Form	Liquid
Buffer	4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human TWIST1/2.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04942 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX04942 Twist1/2 antibody.

Antigen retrieval : Heat mediation with EDTA buffer (pH 8.0).

Dilution : 2 μg/ml


GTX04942 FCM Image


FACS analysis of 4% PFA-fixed Raw264.7 cells using GTX04942 Twist1/2 antibody.

Blue : Primary antibody

Green : Isotype control

Red : Secondary antibody


Dilution : 1 μg/1x10⁶ cells


GTX04942 IHC-P Image

IHC-P analysis of rat brain tissue using GTX04942 Twist1/2 antibody.

Antigen retrieval : Heat mediation with EDTA buffer (pH 8.0).

Dilution : 2 μg/ml


GTX04942 FCM Image

FACS analysis of 4% PFA-fixed RH35 cells using GTX04942 Twist1/2 antibody.

Blue : Primary antibody

Green : Isotype control

Red : Secondary antibody

Dilution : 1 μg/1x10⁶ cells



For full product information, images and publications, please visit our [website](https://www.genetex.com).