

SPNS2 antibody, N-term

Cat. No. GTX04943

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human

Package

400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:100
FCM	1:25

Note : Permeabilization : 90% methanol

Not tested in other applications.

Calculated MW 58 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 68-94 amino acids from the N-terminal region of human SPNS2.
Purification	Purified by protein A column followed by antigen-affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04943 WB Image

WB analysis of various samples using GTX04943 SPNS2 antibody, N-term.

Lane 1 : Human heart tissue lysate


Lane 2 : Human kidney tissue lysate

Lane 3 : Human skeletal muscle tissue lysate

Lane 4 : SK-BR-3 whole cell lysate

Loading : 20µg

Dilution : 1:1000


GTX04943 WB Image

WB analysis of various samples using GTX04943 SPNS2 antibody, N-term.


Lane 1 : Human brain tissue lysate

Lane 2 : 293 whole cell lysate

Lane 3 : SK-BR-3 whole cell lysate

Loading : 20µg

Dilution : 1:1000


GTX04943 WB Image

WB analysis of MDA-MB-231 whole cell lysate using GTX04943 SPNS2 antibody, N-term.

Loading : 35µg



For full product information, images and publications, please visit our [website](http://www.genetex.com).