

GAL3ST1 antibody

Cat. No. GTX04952

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human, Mouse

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
IHC-P	1:50-1:100
FCM	1:10-1:50

Not tested in other applications.

Calculated MW 49 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide between 88-116 amino acids from the Central region of human GAL3ST1.
Purification	Purified by protein A column followed by antigen-affinity chromatography.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX04952 WB Image

WB analysis of mouse cerebellum tissue lysate using GTX04952 GAL3ST1 antibody.
Loading : 35µg


GTX04952 FCM Image

FACS analysis of CEM cells using GTX04952 GAL3ST1 antibody.
Green line : Primary antibody
Blue line : Negative control cell


GTX04952 WB Image

WB analysis of various samples using GTX04952 GAL3ST1 antibody.
Lane 1 : Human kidney tissue lysate
Lane 4 : Mouse kidney tissue lysate
Loading : 20µg
Dilution : 1:2000



For full product information, images and publications, please visit our [website](https://www.genetex.com).