

Chicken Ovalbumin protein

Cat. No. GTX04953-pro

Species Chicken**Package**
1 mg

Properties


Form	Lyophilized powder
Buffer	Reconstitute with 10mM PBS to 0.1-1.0 mg/mL (Do not vortex). Lyophilized from 10mM PBS, 5% trehalose.
Storage	Store at 4°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Region/Sequence	Native Ovalbumin (OVA) protein
Expression System	Egg Whites
Purification	Purified by salt co-precipitation and ionic-Exchange chromatography
Purity	>90% as determined by SDS-PAGE.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

**GTX04953-pro Image**

SDS-PAGE of GTX04953-pro Chicken Ovalbumin protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).