

MDA-MB-361 whole cell lysate

Cat. No. GTX04954

Species Human

Package
100 µg

PRODUCT

Summary

MDA-MB-361 lysate was prepared by homogenization in modified RIPA buffer. Cell debris was removed by centrifugation. Protein concentration was determined with Bio-Rad protein assay. This lysate is for use in Western blotting, 10 µg to 20 µg per lane is recommended for mini gel.

Tissue type

Tumor

Properties

Form Liquid

Buffer 50 mM Tris-HCl pH 6.8, 12.5% glycerol, 1% sodium dodecylsulfate, 0.01% bromophenol blue, 5% β-mercaptoethanol.

Preservative No preservatives

Storage This whole cell lysate is stable for three months at -20°C

Concentration 2 mg/ml (Please refer to the vial label for the specific concentration.)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).