

Cyclin E1 antibody [MSVA-898R] HistoMAX™

Cat. No. GTX04969

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package

500 µl, 100 µl

PRODUCT

Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

IHC-P	1:50 - 1:100
-------	--------------

Note : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Tris-EDTA-based Target Retrieval Solution buffer.

Not tested in other applications.

Product Note

Highly recommended for IHC-P in human tissues.

Properties

Form	Liquid
Buffer	10mM PBS, BSA (Please contact us for sodium azide-free or sodium azide and BSA-free format)
Preservative	sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.2 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human Cyclin E1 protein fragment
Purification	Protein A/G purified
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx04969 IHC-P Image


IHC-P analysis of human placenta (mature) tissue using GTx04969 Cyclin E1 antibody [MSVA-898R] HistoMAX™.

Mature placenta showing a strong Cyclin E1 staining of most cytotrophoblast cells.


GTx04969 IHC-P Image

IHC-P analysis of human ovarian clear cell carcinoma (OCCC) tissue using GTx04969 Cyclin E1 antibody [MSVA-898R] HistoMAX™.

Ovarian clear cell carcinoma with strong Cyclin E1 staining of all tumor cells.


GTx04969 IHC-P Image

IHC-P analysis of human skin tissue using GTx04969 Cyclin E1 antibody [MSVA-898R] HistoMAX™.

Skin with a gradually increasing Cyclin E1 staining of squamous epithelial cells towards the surface cell layers.



For full product information, images and publications, please visit our [website](https://www.genetex.com).