

## NHE3 antibody

Cat. No. GTX04975

|              |            |
|--------------|------------|
| Host         | Rabbit     |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | WB         |
| Reactivity   | Human      |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | Assay dependent      |

Not tested in other applications.

Calculated MW 93 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 2 % Sucrose                                                                                                                                                                                                           |
| Preservative  | 0.09% Sodium Azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 0.5 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                 |
| Immunogen     | A synthetic peptide directed towards the middle region of human SLC9A3:<br>AEDMVTHHTLQYLYKPRQEYKHLYSRHELTPTDEKQDREIFHRTMRKR                                                                                                |
| Purification  | Purified by affinity chromatography                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX04975 WB Image**

WB analysis of RPMI-8226 whole cell lysate using GTX04975 NHE3 antibody.

Dilution : 1 µg/mL


**GTX04975 WB Image**

WB analysis of THP-1 whole cell lysate using GTX04975 NHE3 antibody.

Dilution : 1 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).