

## RHBG antibody

Cat. No. GTX04978

|              |                |
|--------------|----------------|
| Host         | Rabbit         |
| Clonality    | Polyclonal     |
| Isotype      | IgG            |
| Applications | WB, IHC-P, FCM |
| Reactivity   | Human, Mouse   |

Package  
400 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:1000               |
| IHC-P              | 1:50-1:100           |
| FCM                | 1:10-1:50            |

Not tested in other applications.

Calculated MW 47 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS                                                                                                                                                                                                                        |
| Preservative  | 0.09% Sodium Azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | KLH conjugated synthetic peptide between 423-450 amino acids from the C-terminal region of human RHBG.                                                                                                                     |
| Purification  | Purified by protein A followed by antigen-affinity chromatography.                                                                                                                                                         |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX04978 FCM Image**

FACS analysis of K562 cells using GTX04978 RHBG antibody.


Green line : Primary antibody

Blue line : Negative control


**GTX04978 WB Image**

WB analysis of K562 whole cell lysate using GTX04978 RHBG antibody. <

Loading : 35 µg


**GTX04978 IHC-P Image**

IHC-P analysis of mouse kidney tissue using GTX04978 RHBG antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).