

## Mesothelin antibody [YP158] (ready-to-use)

Cat. No. GTX04979

|              |            |
|--------------|------------|
| Host         | Rabbit     |
| Clonality    | Monoclonal |
| Isotype      | IgG        |
| Applications | IHC-P      |
| Reactivity   | Human      |

Package  
6 ml

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| IHC-P              | neat                 |

**Note : Recommendations for antigen retrieval : EDTA Buffer pH 8.0**

Not tested in other applications.

## Properties

|               |                                                                                         |
|---------------|-----------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                  |
| Buffer        | Tissue culture supernatant                                                              |
| Preservative  | 0.09% Sodium Azide                                                                      |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)        |
| Immunogen     | Synthetic peptide corresponding to full-length human mesothelin protein                 |
| Purification  | Unpurified                                                                              |
| Conjugation   | Unconjugated                                                                            |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX04979 IHC-P Image**

IHC-P analysis of human mesothelioma tissue using Mesothelin antibody [YP158] (ready-to-use).  
Formalin fixed paraffin embedded mesothelioma with Mesothelin.



For full product information, images and  
publications, please visit our [website](#).