

Methylated Lysine antibody

Cat. No. GTX04989

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IP, ELISA
Reactivity	Species independent

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:20000
ICC/IF	1:200
IP	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Methylated KLH Conjugated
Purification	Purified by Protein A
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTX04989 WB Image

WB analysis of BSA conjugated with/without Methylated Lysine using GTX04989 Methylated Lysine antibody.

Dilution : 1:1000



GTX04989 ICC/IF Image

ICC/IF analysis of 2% PFA-fixed HeLa cells using GTX04989 Methylated Lysine antibody.

Green : Primary antibody

Blue : DAPI

Dilution : 1:50



For full product information, images and publications, please visit our [website](#).