

SCUBE2 antibody

Cat. No. GTX04992

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:10-1:50
IHC-P	1:50-1:100

Not tested in other applications.

Calculated MW 110 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide between 251-280 amino acids from the Central region of human SCUBE2.
Purification	Purified by protein A followed by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX04992 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX04992 SCUBE2 antibody.




GTX04992 WB Image

WB analysis of various samples using GTX04992 SCUBE2 antibody.

Loading : 20µg

Dilution : 1:1000



GTX04992 ICC/IF Image

ICC/IF analysis of ZR-75-1 cells using GTX04992 SCUBE2 antibody.

Green : Primary antibody

Blue : DAPI



For full product information, images and publications, please visit our [website](https://www.genetex.com).