

ITIH2 antibody

Cat. No. GTX04997

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:50-1:100
FCM	1:10-1:50

Not tested in other applications.

Calculated MW 106 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A KLH conjugated synthetic peptide from the C-terminal region of human ITIH2 (858-886 aa).
Purification	Purified by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

DATA IMAGES

GTx04997 FCM Image


FACS analysis of HepG2 cells using GTx04997 ITIH2 antibody.

Green line : Primary antibody

Blue line : Negative control


GTx04997 IHC-P Image

IHC-P analysis of human hepatocarcinoma tissue using GTx04997 ITIH2 antibody.


GTx04997 WB Image

WB analysis of HepG2 whole cell lysate using GTx04997 ITIH2 antibody.

Loading : 35μg



For full product information, images and publications, please visit our [website](https://www.genetex.com).