

Lactic acid-Lysine antibody

Cat. No. GTX05002

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Species independent

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 50%Glycerol
Preservative	0.02% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a sequence containing Lactic acid-Lysine.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx05002 WB Image

WB analysis of NIH/3T3 cells with/without 100mM Sodium lactate for 24hr treatment using GTx05002 Lactic acid-Lysine antibody.

Loading : 25μg

Dilution : 1:400




GTx05002 WB Image

WB analysis of C6 cells with/without 100mM Sodium lactate for 24hr treatment using GTx05002 Lactic acid-Lysine antibody.

Loading : 25μg

Dilution : 1:400



GTx05002 WB Image

WB analysis of MCF7 cells with/without 100mM Sodium lactate for 24hr treatment using GTx05002 Lactic acid-Lysine antibody.

Loading : 25μg

Dilution : 1:400



For full product information, images and publications, please visit our [website](https://www.genetex.com).