

APC antibody [ARC53932]

Cat. No. GTX05003

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
IHC-P	1:50-1:200

Not tested in other applications.

Calculated MW 312 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA, 50% Glycerol
Preservative	0.05% Proclin300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2600-2760 of human APC (NP_000029.2).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX05003 WB Image

WB analysis of various samples using GTX05003 APC antibody [ARC53932].

Loading : 25µg


Dilution : 1:500


GTX05003 IHC-P Image

IHC-P analysis of mouse colon tissue using GTX05003 APC antibody [ARC53932].

Antigen retrieval : High pressure antigen retrieval performed with 0.01M Citrate Bufferr (pH 6.0)

Dilution : 1:200


GTX05003 IHC-P Image


IHC-P analysis of rat brain tissue using GTX05003 APC antibody [ARC53932].

Red : Primary antibody

Blue : DAPI

Antigen retrieval : Microwave antigen retrieval performed with 0.01M Citrate Buffer (pH 6.0)

Dilution : 1:200


GTX05003 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX05003 APC antibody [ARC53932].

Red : Primary antibody

Blue : DAPI

Antigen retrieval : Microwave antigen retrieval performed with 0.01M Citrate Buffer (pH 6.0)

Dilution : 1:200



For full product information, images and publications, please visit our [website](http://www.genetex.com).