

Mouse Anti-Rabbit IgG antibody [2A9]

Cat. No. GTX05005

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, FCM, ELISA
Reactivity	Rabbit

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
FCM	Assay dependent
ELISA	Assay dependent

Note : For flow cytometry, the suggested use of these reagents is in a final volume of 100µL

Not tested in other applications.

Product Note Minimal reactivity to human, mouse, rat, goat, bovine, horse, donkey, sheep, guinea pig, and chicken serum.

Properties

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant rabbit IgG Fc fragment
Purification	Purified IgG1
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTx05005 ELISA Image

ELISA analysis of serially diluted rabbit IgG and IgM using GTx05005 Mouse Anti-Rabbit IgG antibody [2A9].



For full product information, images and publications, please visit our [website](https://www.genetex.com).