

Trk antibody

Cat. No. GTX05014

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:50-1:200

Note : Antigen retrieval : EDTA Buffer pH 8.0

Not tested in other applications.

Properties

Form	Liquid
Buffer	Batch dependent
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide-KLH conjugate with human Trk A (777-796aa)
Purification	Purified IgG From serum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx05014 IHC-P Image

IHC-P analysis of human cerebellum tissue using GTx05014 Trk antibody.



For full product information, images and publications, please visit our [website](#).