

IRTA1 antibody [A1]

Cat. No. GTX05016

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	IHC-Fr, FCM
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-Fr	Assay dependent
FCM	2-5 µg/ml

Note : Do not need a permeabilization step.

Not tested in other applications.

Product Note This antibody recognizes an epitope within extracellular domain of CD307d.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	DNA-immunization followed by a boost with CD307d-transfected cells
Purification	Purified by Protein A
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

DATA IMAGES

GTX05016 FCM Image

FACS analysis of human lymphocytes using GTX18229 CD19 antibody [LT19] (PE) and GTX05016 IRTA1 antibody [A1].

GTX18229 : 20 µl reagent / 100 µl of peripheral whole blood

GTX05016 : 15 µg/ml (GAM-APC)



For full product information, images and publications, please visit our [website](https://www.genetex.com).