

IRTA1 antibody [A1] (PE)

Cat. No. GTX05016-08

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	FCM
Reactivity	Human

Package
100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	10 μ l reagent in 100 μ l of whole blood or 10 ⁶ cells in a suspension

Note : Do not need a permeabilization step.

Not tested in other applications.

Product Note This antibody recognizes an epitope within extracellular domain of CD307d.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	DNA-immunization followed by a boost with CD307d-transfected cells
Purification	Purified by size-exclusion chromatography
Conjugation	Phycoerythrin (PE) Wavelength


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX05016-08 FCM Image

FACS analysis of human lymphocytes using GTX05016-08 IRTA1 antibody [A1] (PE) and GTX18224 CD19 antibody [LT19] (APC).

GTX05016-08 : 10 μ l reagent in 100 μ l of peripheral whole blood

GTX18224 : 4 μ l reagent in 100 μ l of peripheral whole blood



For full product information, images and publications, please visit our [website](https://www.genetex.com).