

SLC25A6 antibody

Cat. No. GTX05019

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human

Package
400 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:10-1:50
IHC-P	1:10-1:50
FCM	1:10-1:50

Not tested in other applications.

Calculated MW 33 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide from the Central region of human SLC25A6(128-155 aa).
Purification	Purified by protein A followed by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX05019 FCM Image

FACS analysis of A549 cells using GTX05019 SLC25A6 antibody.


Green line : Primary antibody

Blue line : Negative control



GTX05019 IHC-P Image

IHC-P analysis of human liver tissue using GTX05019 SLC25A6 antibody.



GTX05019 ICC/IF Image

ICC/IF analysis of A549 cells using GTX05019 SLC25A6 antibody.

Green : Primary antibody

Blue : DAPI



For full product information, images and publications, please visit our [website](#).