

PRR11 antibody [OTI3C8]

Cat. No. GTX05021

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:2000
IHC-P	1:150

Not tested in other applications.

Calculated MW 40 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 1% BSA, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein fragment (125-360 aa) corresponding to human PRR11 (NP_060774) produced in E.coli.
Purification	Purified by Protein A/G
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX05021 IHC-P Image

IHC-P analysis of human thyroid tissue using GTX05021 PRR11 antibody [OTI3C8].
Antigen retrieval : Heat mediated antigen retrieval step in 10mM citric buffer (pH6.0) at 100°C




GTX05021 IHC-P Image

IHC-P analysis of human thyroid tissue using GTX05021 PRR11 antibody [OTI3C8].
Antigen retrieval : Heat mediated antigen retrieval step in 10mM citric buffer (pH6.0) at 100°C



GTX05021 WB Image

WB analysis of various samples using GTX05021 PRR11 antibody [OTI3C8].
Loading : 35µg



GTX05021 WB Image

WB analysis of HEK293T cells with/without PRR11 overexpression using GTX05021 PRR11 antibody [OTI3C8].
Loading : 5µg



For full product information, images and publications, please visit our [website](https://www.genetex.com).