

Goat Anti-Human IgG antibody, pre-adsorbed (FITC)

Cat. No. GTX05041

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IHC-Fr, FCM, IP, ELISA, Depletion, Multiplexing
Reactivity	Human

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	$\leq 1 \mu\text{g}/10^6$ cells
IP	Assay dependent
ELISA	1:100-1:400
Depletion	Assay dependent
Multiplexing	Assay dependent

Not tested in other applications.

Product Note Preadsorbed with human IgM and IgA. May react with immunoglobulins from other species.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified by affinity chromatography From serum
Conjugation	Fluorescein isothiocyanate (FITC) Wavelength




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX05041 ELISA Image

ELISA analysis of purified human IgG, IgM and IgA using serially diluted GTX05041 Goat Anti-Human IgG antibody, pre-adsorbed (FITC).



For full product information, images and publications, please visit our [website](https://www.genetex.com).