

Pan phosphoserine/threonine antibody [ARC64485]

Cat. No. GTX05060

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Species independent

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
IHC-P	1:50-1:200

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA, 50% Glycerol
Preservative	0.09% Sodium Azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a sequence containing phosphorylated Serine/Threonine.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTK05060 IHC-P Image

IHC-P analysis of human thyroid tissue using GTK05060 Pan phosphoserine/threonine antibody [ARC64485].

Antigen retrieval : High pressure mediated antigen retrieval performed with 0.01 M citrate buffer (pH 6.0)


Dilution : 1:200


GTK05060 IHC-P Image

IHC-P analysis of mouse testis tissue using GTK05060 Pan phosphoserine/threonine antibody [ARC64485].

Antigen retrieval : High pressure mediated antigen retrieval performed with 0.01 M citrate buffer (pH 6.0)

Dilution : 1:200


GTK05060 IHC-P Image

IHC-P analysis of mouse colon tissue using GTK05060 Pan phosphoserine/threonine antibody [ARC64485].

Antigen retrieval : High pressure mediated antigen retrieval performed with 0.01 M citrate buffer (pH 6.0)

Dilution : 1:200


GTK05060 IHC-P Image

IHC-P analysis of human colon tissue using GTK05060 Pan phosphoserine/threonine antibody [ARC64485].

Antigen retrieval : High pressure mediated antigen retrieval performed with 0.01 M citrate buffer (pH 6.0)

Dilution : 1:200



For full product information, images and publications, please visit our [website](https://www.genetex.com).